



Campylobacter Growth Supplement (Liquid) SR0232E



10 x :	Sodium pyruvate	0.125 g	Sodium metabisulphite	0.125 g
	Ferrous sulphate	0.125 g		

Precautions: Use sterile techniques at all times. Do not use beyond stated expiry date.
 Precauciones: Manipular asépticamente durante todo el procedimiento. No utilizar una vez caducado.
 Vorsichtsmassnahmen: Nur unter aseptischen Bedingungen arbeiten. Nicht nach Ablauf des Verfallsdatum verwenden.
 Précaution: Travailler à chaque fois de manière stérile. Ne pas utiliser après la date d'expiration.
 Precauzioni: Operare sempre in condizioni di sterilità. Non utilizzare dopo la data di scadenza.
 Forholdsregel: Anvend steril teknik til hver en tid. Må ikke anvendes ud over den angivne udløbsdato.
 Försiktighet: Använd alltid sterila tekniker. Använd inte produkterna efter deras utgångsdatum.
 Precauções: Utilizar sempre a técnica estéril. Não utilizar depois de ter caducado o prazo de validade indicado.
 Προφυλάξεις: Να χρησιμοποιείτε πάντοτε άσηπτη τεχνική. Μην χρησιμοποιείτε το προϊόν μετά την αναγραφόμενη ημερομηνία λήξης.
 注意: 常に滅菌処理を行ってください。記載されている使用期限を過ぎたものは、絶対に使用しないでください。



500 ml 50°C

Sterile medium	Sterilit medie
Del medio estéril	Sterilit medium
Steriles Medium	Meio estéril
Milieu stérile	Στείρο μέσο
Di coltura sterile	滅菌した培地



Campylobacter Agar Base (CM0689)
 or Columbia Agar Base (CM0331)
 or Blood Agar Base No. 2 (CM0271)
 or Nutrient Broth No. 2 (CM0067)

+

5-7% v/v lysed horse blood (SR0048)
 or 5-7% v/v defibrinated horse blood (SR0050)
 or 5-7% v/v defibrinated sheep blood (SR0051)

+

Campylobacter Selective Supplement
 (SR0069 or SR0085 or SR0098 or SR0117
 or SR0204 or SR0214)



Mix well
 Mezclar bien
 Gut mischen
 Bien mélanger
 Mescolare con cura

Bland godt
 Blanda ordentligt
 Misturar bien
 Αναμείξτε καλά
 良く混和します



Dispense into sterile containers
 Dispensar en envases estériles
 In sterile Behältnisse abfüllen
 Distribuire in contenitori sterili
 Répartir de façon stérile

Dispensér i sterile beholdere
 Dispensera i sterila behållare
 Dispensar em recipientes estéreis
 Μεταφέρετε σε στείρους περιέκτες
 滅菌済みの容器に分注してください